

# SPECIFICATION

## One-pc open Pouch (Pediatrics velcro closure)



20 pouches, per box

### Product Feature:

1. Pouch size: 13.5X23cm
2. Pouch volume: 350ml
3. Using time: 5-7 days for colostomy,  
3-5 days for ileostomy,
4. Pouch film: 5 layer co-extruded EVA film,  
odor-free and noiseless  
thickness: 0.075-0.085mm  
transparent/beige
5. Skin barrier: Diameter=9x9cm,  
Thickness: 0.75mm  
Cutting range: 15-50mm
6. With bar filter
7. Single side nonwoven liner

### Product Safety :

The permanence glutinosity test:	Pass
Skin irritation test:	Pass
Cytotoxicity test:	Pass
Delayed hypersensitivity test:	Pass
Anti-leakage test:	Pass
Odor impermeable test:	Pass



### SPECIAL STORAGE CONDITIONS:

Store in a cool dry place, in which air containing high levels of oxygen/ozone away from direct heat and sunlight.

**SHELF LIFE: 3 years if stored as recommended above**

### Company information:

BAO-Health Medical Instrument Co., Ltd.  
2nd Floor of 8-2 Building, Gaoqiao Industrial Zone,  
Tongxiang, Zhejiang, P.R.China, 314515

T:+86-573-80889033 F:+86-573-80889034  
www.bao-health.com info@bao-health.com

